


APPLICATION NOTE

Multi-parametric phenotyping of *C. elegans* neuronal degeneration

PURPOSE

Neuro-degeneration studies in *C. elegans* require high-resolution imaging of a large number of immobilized animals. The vivoChip® enables detection of sub-lethal, multi-parametric, and neuron-specific phenotypic changes in 40 *dat-1::gfp* labeled animals at once.



METHOD

- Up to 40 adult animals are immobilized in the vivoChip®-2x within 3 min.
- Dopaminergic neurons labeled with *dat-1::gfp* are imaged at high magnifications (20x, 0.75NA).
- Neuronal degenerations are induced in adult animals by treatment with the neurotoxin MPP+ iodide.
- Defined classes of structural aberrations (phenotypes) are evaluated in the immobilized animals.




Figure 1. Images of BY200 (*dat-1::gfp*)¹ animals immobilized in the vivoChip®-2x.

| Images from vivoChip®-2x | ImageJ analysis | Subtle phenotypes |
|--------------------------|-----------------|---|
| | | <ul style="list-style-type: none"> • Beading (BD) • Breaks (BR) |
| | | <ul style="list-style-type: none"> • Beading (BD) • Tip loss (TL) |
| | | <ul style="list-style-type: none"> • Beading (BD) • Breaks (BR) • Deformation (DF) |
| | | <ul style="list-style-type: none"> • Beading (BD) • Breaks (BR) |

CONCLUSIONS

- The vivoChip®-2x enables rapid immobilization and high-resolution imaging of entire neuronal processes using air objectives (20x, 0.75NA) or oil objectives (up to 100x).
- Simultaneous immobilization of 40 animals side-by-side facilitates fast imaging and simple manual scoring of multi-parametric neuronal degeneration phenotypes using ImageJ software.
- High-resolution imaging of anesthetized animals (for complete immobilization) allows researchers to identify even subtle neuronal degeneration phenotypes enabling investigation of the mechanisms of neuronal toxicity.

¹BY200 (*dat-1::gfp*) strain was provided by Dr. Michael Aschner, Albert Einstein College of Medicine.